Recombinant Human Transferrin

Optiferrin™

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
Stem Cell & Regenerative Medicine Cell Lines

- Mesenchymal
- Embryonic
- iPS
- Cardiac
- Liver
- Kidney
- Keratinocyte
- Fibroblast
- Red blood
- Neural
- Osteoblast
- T cells

Research and Production Cell Lines

- CHO
- Hybridoma
- Vero
- HEK293
- Rat crypt
- HT29
- MDCK
- Hamster kidney
- PER.C6

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




CD34+ stem cells (HSC) isolated from human umbilical-cord mononuclear cells blood were expanded in culture for 6 days. In basal media and 1 g/L CellSiam S Recombinant Human Albumin, different iron sources and iron carriers were supplemented at a concentration of 50 mg/L. The relative effect of supplementation on proliferation compared to conditions without the inclusion of any iron carrier or other iron source can be observed. Cell count was normalized using Plasma Human Serum Albumin and Iron Chelate (Plasma HSA/Iron Chelate = Normalized to 100).

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Functional lentivirus titer is increased with excess transferrin. HEK-293T cells were transfected with 2nd generation lentivirus plasmids and total lentivirus vectors were harvested 48 hours later and titered by functional transduction assay with HT-1080 cells analyzed by flow cytometry. Error bars represent standard deviation. OptiPEAK HEK293t is abbreviated as OptiHEK, transferrin is abbreviated as Tf. Estimated Infectious Units (IFU)/ mL is the average functional titer from three individual experiments with two technical replicates each.

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Transfection efficiency of HEK-293T is improved with excess transferrin. HEK-293T cells were transfected with a GFP reporter plasmid and analyzed 48 hours later by flow cytometry. Error bars represent standard deviation, * and ** indicate statistical significance. OptiPEAK HEK293t is abbreviated as OptiHEK. Mean percent GFP positive cells was calculated from three individual experiments with two technical replicates each.

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J Inorg Biochem. 2012 Nov;116:37-44. doi: 10.1016/j.jinorgbio.2012.07.005. Epub 2012 Jul 11. Biochemical and structural characterization of recombinant human serum transferrin from rice (*Oryza sativa* L).
Protein Expr Purif. 2010 Nov;74(1):89-79. doi: 10.1016/j.pep.2010.04.019. Epub 2010 May 4. Expression, purification, and characterization of recombinant human transferrin from rice (*Oryza sativa* L).

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