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New Human Derived Product Advanced for Cell Culture and Stem Cell therapy

AventaCell BioMedical Co., Ltd. ("AventaCell") is among the world leaders devoted to developing novel human-derived products for use in cell culture and tissue regeneration. Helios Bioscience is AventaCell's product family brand for their products to be used in cell culture and tissue regeneration. AventaCell technologies used in the Helios Bioscience line offer new human-derived solutions for use in translational research of cell and tissue-based therapies to meet the need for animal serum-free cell expansion and production. The demand for safe, efficient and cost-effective cell expansion and production is rapidly increasing with the growth in cell therapy and regenerative medicine research and clinical development. Helios Bioscience products are designed to support expansion and production of a broad range of cells including mesenchymal stem cells and multiple immunocell lines. AventaCell is committed to providing animal serum-free products to accelerate the research, development and commercialization of safe and efficacious cell and tissue-based therapeutics.

Non-Xenogenic Supplement
Animal Serum Free

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
Figure 1. T cell expanded from PBMC in medium (Culture media supplier B) supplemented with UG-P GI or autologous plasma. (A) Total cell fold expansion (B) Percentage of live cells.

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
Figure 2. NK cell expanded from PBMC in medium (culture media supplier C) supplemented with UG-P GI or AB serum. (A) Total cell fold expansion. (B) Cytotoxicity of NK cells against K562 cells at 3:1, 5:1 effectors : target ratios. (C) Percentage of live cells, Percentage of NK cells, Specific expansion fold of NK cells.

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
Figure 3. Expansion of NK cell line (NK-92) in various cell culture media with IL-2 (200U/mL) and different supplements. After 12 days in culture, expanded cells were analyzed for (A) Expansion fold of NK cells (B) Cytotoxicity of NK cells against K562 cells at 5:1 effectors : target ratio.

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